

STERILIZATION OF OBJECTS USING H₂O₂

Bob Dotterer, VP of Engineering; Alan Campbell, VP of Marketing

Abstract

Caron's patented "dry" vaporized hydrogen peroxide system can be validated to sterilize large volumes of even thermally sensitive items that aren't compatible with traditional high heat sterilization.

Introduction

The recent COVID-19 crisis has highlighted the need to be able to decontaminate mass quantities of a wide range of reusable user-contact products, such as face shields, masks, and other personal protective equipment, for emergency use. Most of these products incorporate synthetic materials that are not suitable for repeated heat-based decontamination or sterilization. There is also a need to decontaminate or sterilize other types of heat-sensitive employee or user-contact products, such as filled sample vials. While low temperature decontamination technologies such as Ethylene Oxide (EtO), Peracetic Acid and Paraformaldehyde exist to fill this need, they have long-term toxicity issues that have limited their use, mostly to specialized offsite facilities. Hydrogen Peroxide (H₂O₂) is another low-temperature option. This powerful oxidizer has well documented and accepted anti-microbial properties, permitting a true Sterilization claim. Since it is fast-acting and catalyzes quickly, it offers short cycle times. Without the need to heat a chamber up to a high temperature, it also minimizes stress on unit components, sensors, and electronics. Unlike previous generations of low temperature gas sterilants, H₂O₂ decomposes into oxygen and water vapor after use, producing no lingering toxic chemicals or byproducts. As a result, H₂O₂ has become a major contamination elimination technology in both clinical sterile processing departments and life science research.

Caron has advanced the application of vaporized H₂O₂ sterilization technology by making it much simpler to use, and available within a large volume enclosure. Caron's gFog™ H₂O₂ technology pairs reduced user involvement with the shortest cycle duration. Caron integrates a dehumidification process into the cycle conditioning phase (US patent number 10,738,271). Unlike simple "wet" H₂O₂ vapor cycles, which saturate the chamber with an H₂O₂ mist, creating sterilant condensation and resulting sticky residue, Caron's dry process injects hydrogen peroxide vapor into the airstream in a controlled method, keeping the

amount of H₂O₂ and water vapor below saturation point, usually around 90% RH. A humidity sensor is used to monitor the amount of H₂O₂ in the air, and electronic control systems throttle vapor injection and maintain it at a constant level. By keeping the humidity below saturation point, no condensation forms anywhere within the chamber. Instead of a process that requires interior disassembly, drying, and reassembly after use, Caron's dry vapor cycle process permits immediate unit use post-cycle, minimizing both downtime and effort.

Background

Caron's cycle was designed to sterilize non-porous chamber surfaces, with sufficient capacity and duration to also decontaminate or sterilize high volumes of customer-provided porous or non-porous product within the chamber.

Caron's H₂O₂ process consists of three phases:

- Temperature (dehumidification)/ H₂O₂ increase to necessary levels (conditioning)
- Hold Temperature & H₂O₂ levels (kill phase)
- H₂O₂ removal (inactivation)

Sterilant volume and cycle duration parameters were determined through both theoretical and empirical test data:

- A 35% H₂O₂ concentration was selected for two reasons:
 - A higher concentration sterilant solution decreases the total volume of liquid that must be injected into the chamber and reduces potential for condensation.
 - It is the highest commonly used & commercially available concentration that is still under the US Department of Transportation threshold for passenger, cargo air freight, and rail shipment.
- H₂O₂ exposure D-value for chamber surfaces was established during an initial round of development and testing. Refer to "Sterilization of Reach-In Incubators Using H₂O₂" for specific details.
- Injecting H₂O₂ during the conditioning phase is critical to achieving a rapid sterilization cycle. 20% of the sterilization exposure is achieved during the conditioning phase, shortening the sterilization phase. (patented)
- Additionally, the inactivation phase begins with an existing high level of H₂O₂ concentration in the chamber. As this H₂O₂ level decreases over the 60 min time duration to a safe level, the H₂O₂ continues to act as a sterilant and provide further margin of product exposure.
- *Bacillus stearothermophilus* is the most prevalent Biological Indicator (BI) organism for validating H₂O₂ because it is one of the most resistant to hydrogen peroxide.

Caron assumes a unit temperature start value of 25°C prior to cycle start:

- If the chamber has been 'off' for hours, the internal temperature will be ambient, typically 18°C-22°C. In this case, 25°C is also a conservative cycle start value number.
- If the chamber had run through a cycle recently, the walls may be above 25°C, but the air temperature will be near room temperature, because the doors must be opened to load customer product and initiate the test.

Prior testing with BI's proved that a chamber cycle could achieve a complete 12-log reduction of all BI's by ramping from 25°C back up to 37°C and holding that temperature for 18 minutes while injecting H₂O₂ to 90% RH. Caron's Decontamination Chamber features a 45°C cycle, to further reduce cycle duration.

Caron's H₂O₂ injection takes place via a durable heated probe, which penetrates a disposable container, contacting the liquid sterilant within. The probe vaporizes H₂O₂, which is then drawn into the blower fan and circulated through the chamber via a pressurized back duct and uniform horizontal airflow across each shelf. H₂O₂ injection can only be initiated through the designated process, which includes an electronic door lock function, eliminating the potential for H₂O₂ release outside of the sealed chamber environment.

Caron uses a large volume activated charcoal filter to catalyze H₂O₂ into H₂O and O₂. This long-life filter is designed for continuous unit operation. It has been determined that a process inactivation step of 60 minutes will lower the H₂O₂ to a safe level.

Materials and Methods


Given the varied materials, quantities, and sample construction techniques possible for products being sterilized, standard BI's may not be appropriate. Porous materials, such as textiles or paper, and samples incorporating narrow crevices or interior surfaces will require longer exposure and greater sterilant volume within the cycle to produce the desired result. It is possible that certain material volumes or configurations may not be sterilizable within the standard unit cycle. Consult with a recognized third-party test laboratory, such as Mesa Laboratories (Bozeman, MT), to determine whether a standard test coupon is available for your application, or if a material or configuration specific to your application should be used. Typically, the latter option would involve sending a defined sample to the test lab for inoculation, after which validation would proceed in a similar way to that of standard BI's.

- * 7800-25, Decontamination Chamber
- * STER309 or STER310, 35% Hydrogen Peroxide, container of 56ml (10x or 50x)
- * Either standard Biological Indicator (BI) (reference Apex Discs, Geobacillus stearothermophilus spores, Catalog number APEX-456, by Mesa Labs or samples of target product load type(s), as determined in consultation with a 3rd party test laboratory.
- * Standard Chemical Indicator (CI) strips, 3M Comply Hydrogen Peroxide, Catalog number 1248
- * Data logger, Keysight model 34972A LXI
- * Basic refrigerator, maintain temperature 2-8°C

BI Test Locations Within the Chamber



Summary of H₂O₂ Cycle Steps



* Pre-validated for chamber surfaces, owner validation required for product load.

| Stage | Time | Temperature | Humidity (H ₂ O ₂ + H ₂ O) |
|-----------------------|------------|--------------|---|
| Conditioning | 20 minutes | 25°C to 45°C | Ambient to 90% |
| Sterilization* | 25 minutes | 45°C | 90% |
| Inactivation | 60 minutes | 46°C** | 90% to ambient |

* Pre-validated for chamber surfaces, owner validation required for product load

** Novel inactivation process patent pending

Validation Test Procedure

1. Remove all contents from chamber.
2. Open unit door.
3. Position BI's in representative locations within the chamber, typically back-left of the top shelf, front-right of the bottom shelf, and left center and right center of the middle shelf.
4. Attach a control BI to the exterior wall of the chamber.
5. Attach a CI to the clip located on the right-hand side of the integral H₂O₂ module.
6. Install 35% H₂O₂ canister within the integral H₂O₂ module, and close module door.
7. Close the unit door.
8. Initiate a cycle from the unit touchscreen controller.
9. Verify unit is performing properly throughout the sterilization cycle.
10. When the cycle is complete, confirm H₂O₂ exposure by examining the CI, and then remove all BI's & store in refrigerator at 5°C.
11. Send all exposed BI's to third-party test facility (in this case, MesaLab) for overnight processing.

The no-touch H₂O₂ canister does not require measuring or pouring and can be disposed of after use. Estimated cumulative user time to initiate and conclude a cycle is less than 5 minutes.

BI Processing

12. Per protocol and third-party test lab's internal procedures

Validation Results for Growth

Negative growth results will indicate successful inactivation of that test coupon. The outside control coupon, however, should test positive for growth. If the test coupon fails to test positive, re-run the validation process until a positive result is achieved.

Conclusions

Caron's gFog H₂O₂ cycle achieves fast, documented, and highly reproducible low temperature chemical decontamination for large reach-in chambers. This cycle was developed using established and recognized methods and test materials and employs Caron's uniform directed airflow to ensure consistent heat, humidity, and H₂O₂ distribution throughout the chamber.

With the proper validation testing, consistent sterilization of chamber loads can also be established, even in large quantities, and with porous materials. While this might not apply to every product type, vaporized H₂O₂'s ability to rapidly sterilize surfaces and materials make it an ideal choice for many sample types and applications.

Bibliography and References

B. Dotterer and A. Campbell

Sterilization of Reach-In Incubators Using H₂O₂

Kokubo, M., T. Inoue, and J. Akers.

Resistance of common environmental spores of the genus *Bacillus* to vapor hydrogen peroxide vapor.

PDA J. Pharm. Sci. Technol. 52: 228-231, 1998.

ANSI/AAMI/ISO 11138-1

Sterilization of health care products – Biological

Indicators-Part 1: General requirements

Association for the Advancement of Medical Instrumentation, 2006.

S. Radl, S. Ortner, R. Sungkorn, and J. Khinast.

The Engineering of Hydrogen Peroxide Decontamination Systems

J Pharm Innov, 4: 51-62, 2009.

C. Hultman, A. Hill, and G. McDonnell.

The Physical Chemistry of Decontamination with Gaseous Hydrogen Peroxide

Pharmaceutical Engineering, Jan/Feb: 22-32, 2007.

Vaisala Oyj.

Humidity Conversion Formulas

B210973EN-F, 2013.

B. Unger-Bimczok, V. Kottke, C. Hertel.

The Influence of Humidity, Hydrogen Peroxide Concentration, and Condensation on the Inactivation of *Geobacillus*

stearothermophilus Spores with Hydrogen Peroxide Vapor

J Pharm Innov, 3:123-133, 2008.

D. Watling, C. Ryle, M. Parks, and M. Christopher.
Theoretical Analysis of the Condensation of Hydrogen peroxide
Gas and Water Vapour as used in Surface Decontamination
Bioquell Pharma, 2006.

R. Jones, J. Drake, and D. Eagleson
Using H₂O₂ to decontaminate biological safety cabinets
The Baker Company, Acumen Volume 1, No 1, 1993

J. Agalloco, J Akers.
Overcoming Limitations of Vaporized Hydrogen Peroxide
Pharmaceutical Technology, Volume 37, Issue 9, 2013.

J. Spiegelman and D. Alvarez.
Cheating Raoult's Law to Enable Delivery of Hydrogen
Peroxide as a Stable Vapor
Gases & Instrumentation: 14-19, Jan/Feb 2015.

US Department of Transportation.
49CFR 172.101, Hazardous Materials Table
US Code of Federal Regulations, amended 12/21, 1990.