

CO₂/O₂ INCUBATOR / OASIS

Models: 6400 / 6404

Model	6400-1	6400-4	6404-1	6404-4
Temperature Range	5°C above ambient to 60°C			
Temperature Control	±0.1°C at 37°C			
Temperature Uniformity	±0.3°C at 37°C			
Temperature Sensor	Precision Thermistor			
Humidity Range	Elevated up to 95% @ 37°C			
CO ₂ Range	0-20% CO ₂			
CO ₂ Control	±0.1% CO ₂			
CO ₂ Sensor	Infrared CO ₂ Sensor			
O ₂ Range	N/A		1 to 21% O ₂	
O ₂ Control			±0.1% O ₂	
O ₂ Sensor			Fuel Cell	
Interior Dimensions	20"W x 21"D x 25"H (50.8 cm x 53.3 cm x 63.5 cm)			
Interior Construction	Polished Stainless Steel			
Exterior Dimensions	26"W x 26"D* x 36"H (66 cm x 66 cm x 91.4 cm)			
Exterior Construction	Powder Coated Cold Rolled Steel			
Work Space	6 cu. ft. (170 liters)			
# of Shelves	4 Standard; 13 Maximum			
Shelf Construction	Perforated Stainless Steel			
Shelf Dimensions	18.5" x 18.3" (47 cm x 46.5 cm)			
Electrical	115V, 15A, 60 Hz	230V, 8A, 50/60 Hz	115V, 15A, 60 Hz	230V, 8A, 50/60 Hz
Shipping Weight	264 lbs. (120 kg.)			

Specifications are subject to change without notice. *Add 2.75" for handle.