

CO₂ INCUBATOR / SHAKERS + CELL ROLLERS

Models: 7401-25 / 7402-25 / 7403-25

Models: 7401-33 / 7402-33 / 7403-33

Model	7401-25-1 7401-25-2 7401-25-3	7401-33-1 7401-33-2 7401-33-3	7402-25-1 7402-25-2 7402-25-3	7402-33-1 7402-33-2 7402-33-3	7403-25-1 7403-25-2 7403-25-3	7403-33-1 7403-33-2 7403-33-3
Description	CO₂ Incubator for Shakers		CO₂ Incubator for Cell Roller		CO₂ Incubator for Shaker & Cell Roller	
# of Shaker-Ready Shelves	2	3	N/A		1	
Cell Roller Height Capacity	N/A		5-Tiers	7-Tiers	3-Tiers	5-Tiers
Standard Shelf Clearance**	23.5" (59.7 cm)	19" (48.3 cm)	N/A		17.5" (44.5 cm)	
Temperature Range	20°C-60°C (at 200W heat load max.)					
Temperature Sensor	3-wire RTD					
Temperature Control	±0.1°C					
Temperature Uniformity	±0.3°C					
CO ₂ Range	0-20% CO ₂					
CO ₂ Control	±0.1% CO ₂					
CO ₂ Sensor	Infrared CO ₂ Sensor					
Interior Dimensions	32"W x 27"D x 52.7"H (81.3 cm x 68.6 cm x 133.9 cm)	32"W x 27"D x 65.7"H (81.3 cm x 68.6 cm x 166.9 cm)	32"W x 27"D x 52.7"H (81.3 cm x 68.6 cm x 133.9 cm)	32"W x 27"D x 65.7"H (81.3 cm x 68.6 cm x 166.9 cm)	32"W x 27"D x 52.7"H (81.3 cm x 68.6 cm x 133.9 cm)	32"W x 27"D x 65.7"H (81.3 cm x 68.6 cm x 166.9 cm)
Interior Construction	Stainless Steel					
Exterior Dimensions	35.5"W x 33.3"D* x 77.1"H (90.2 cm x 84.6 cm x 195.8 cm)	35.5"W x 33.3"D* x 90.1"H (90.2 cm x 84.6 cm x 228.9 cm)	35.5"W x 33.3"D* x 77.1"H (90.2 cm x 84.6 cm x 195.8 cm)	35.5"W x 33.3"D* x 90.1"H (90.2 cm x 84.6 cm x 228.9 cm)	35.5"W x 33.3"D* x 77.1"H (90.2 cm x 84.6 cm x 195.8 cm)	35.5"W x 33.3"D* x 90.1"H (90.2 cm x 84.6 cm x 228.9 cm)
Exterior Construction	Cold Rolled Steel, Powder Coated					
Shelf Construction	Type 304 Stainless Steel, Electro polished		N/A		Type 304 Stainless Steel, Electro polished	
Shelf Dimensions	29.25"W x 24.45"D (74.3 cm x 62.1 cm)		N/A		29.25"W x 24.45"D (74.3 cm x 62.1 cm)	
Work Space	25 ft ³ (708 Liters)	33 ft ³ (934 Liters)	25 ft ³ (708 Liters)	33 ft ³ (934 Liters)	25 ft ³ (708 Liters)	33 ft ³ (934 Liters)
Specifications are subject to change without notice; are based on 20°C ambient & standard voltage. Shaker and Cell Rollers not included. *Add 2.75" for handle. **Shaker Shelves are set in specific locations standard at the factory, but are designed to easily relocate.						