

HEATED / HUMIDIFIED INCUBATOR

Models: 7002-10 / 7002-25 / 7002-33

Models: 7003-10 / 7003-25 / 7003-33

Model	7002-10 (-1,-2, -3)	7003-10 (-1,-2, -3)	7002-25 (-1,-2, -3)	7003-25 (-1,-2, -3)	7002-33 (-1,-2, -3)	7003-33 (-1,-2, -3)
Temperature Range	10°C above ambient to 70°C					
Temperature Control	±0.1°C					
Temperature Uniformity	±0.3°C					
Humidity Range	N/A	Ambient to 95% RH	N/A	Ambient to 95% RH	N/A	Ambient to 95% RH
Humidity Control	N/A	±3% RH	N/A	±3% RH	N/A	±3% RH
Humidity Sensor	N/A	Capacitive	N/A	Capacitive	N/A	Capacitive
Interior Dimensions	23" W x 25.8" D x 29.8"H (58.4 cm x 65.5 cm x 75.7 cm)		32" W x 27" D x 52.7"H (81.3 cm x 68.6 cm x 133.9 cm)		32" W x 27" D x 65.7"H (81.3 cm x 68.6 cm x 166.9 cm)	
Interior Construction	Type 304, 2B Finish, Solid Stainless Steel					
Exterior Dimensions	45.5" W x 31.7" D* x 37" H (112.3 cm x 82.8 cm x 92.7 cm)		35.5" W x 33.3" D* x 77.1" H (90.2 cm x 84.6 cm x 195.8 cm)		35.5" W x 33.3" D* x 90.1" H (90.2 cm x 84.6 cm x 228.9 cm)	
Exterior Construction	Cold Rolled Steel, Powder Coated					
Work Space	10 cu. ft. (283 liters)		25 cu. ft. (708 liters)		33 cu. ft. (934 liters)	
# of Shelves	3 Standard; 15 Maximum		4 Standard; 25 Maximum		5 Standard; 31 Maximum	
Shelf Construction	Type 304, Perforated Stainless Steel, Electropolished					
Shelf Dimensions	22.5" W x 25.4" D (57.2 cm x 64.5 cm)		29.25" W x 24.45" D (74.3 cm x 62.1 cm)		29.25" W x 24.45" D (74.3 cm x 62.1 cm)	

*Add 2.75" for handle.

Model	7002-10-1	7002-10-2	7002-10-3	7002-25-1	7002-25-2	7002-25-3	7002-33-1	7002-33-2	7002-33-3
	7003-10-1	7003-10-2	7003-10-3	7003-25-1	7003-25-2	7003-25-3	7003-33-1	7003-33-2	7003-33-3
Electrical	115V 60 Hz 12A	208/230V 60 Hz 6A	230V 50 Hz 6A	115V 60 Hz 12A	208/230V 60 Hz 6A	230V 50 Hz 6A	115V 60 Hz 12A	208/230V 60 Hz 6A	230V 50 Hz 6A
Shipping Weight lbs.	420		540**	675		815**	765		915**

*Add 2.75" for handle. **Includes Export Shipping Crate.

NOTE: Specifications are subject to change without notice. Specifications are based on 20°C ambient and standard voltage.